

Hand-operated punch for DBS (Dried Blood Spot) testing

BARINASU



▲ Desk-top Punch

▲ Handheld Punch

▲ Showing how to punch DBS sample by desk-top punch


Punch DBS samples without burrs or contamination

"BARINASU", a hand-operated punch for DBS testing, is used to punch DBS samples for the newborn screening and the genetic screening. The precision of DBS analysis shall be improved because it punches the samples without burrs or contamination.

Features

- Check a punching position visually and punch it
- Easy cleaning with a cotton swab after use
- Hard to get an accretion like blood on the cutting surface
- Light weight and portable
- Easy disassembly and cleaning

Comparison of punched DBS samples (1/8 inch)



Other brand

Nogami brand

Specifications

| | Desk-top Punch | Handheld Punch |
|-----------------|--------------------------------|--------------------------------|
| Dimensions (mm) | 60.0×155.0×109.0 W x L x H | 18.0×117.0×77.0 W x L x H |
| Stamping Size | Cicle(Diameter) 1.00~5.00mm | Cicle(Diameter) 1.00~5.00mm |
| Weight | 1.7kg | 0.17kg |

Stamping Size Punching size is specified by customer

Please talk with us about your problems: such as quality of punched materials.